



Spectrometer-Based Micro-Plate Reader (ELISA) (BMG LABTECH, Germany)

The innovative absorbance microplate reader has the flexibility to perform assays quickly and easily in both microplates or via the built-in cuvette port:

- Captures a full UV/visible spectrum in less than 1 sec/well.
- Quantification of DNA and RNA using small volume 2 μ L using the Lvis plate feature (**Nano-drop like feature**).
- Simple push button operation, and capacity to store individual assay protocols.
- Possible incubation at different temperatures with shaking.
- Comprehensive, and flexible software for analysis of results.

Other available Instruments:

- Rotary evaporator (Buchi, Switzerland)
- Ultrasonic water bath (Elmasonic, Germany)
- Shaking water bath (Thermo Scientific, USA)
- pH meter (Jenway, USA)
- Spinning mini centrifuge for PCR tubes/strips (Thermo Scientific, USA)
- Hot plate magnetic stirrer (STURAT- Cole Palmer, USA)
- Vortex mixer (Thermo Scientific, USA)
- Electric digital balance (Boeco, Germany)
- Bench fee services are available.

Contact us



01001996700



central.lab@pharma.asu.edu.eg



<https://pharma.asu.edu.eg/page/65>



Faculty of Pharmacy – Ain Shams University, Street of African Union Organization, Abbaseya, Cairo, 11566, Egypt.

Chairman of Board of Directors

Prof. Amany Osama Kamel

Director

Associate Prof. Dina Hamada Kassem

Deputy Director

Dr. Hend Zaki Yamani



Central Research Lab

Faculty of Pharmacy, Ain Shams University





Cooling Centrifuge (Eppendorf, Germany)

Multipurpose refrigerated centrifuge with aerosol-tight rotor(s)

Available rotors:

- **Rotor for 1.5/2 ml tubes:**
Max speed: 20,000 × g (14,000 rpm)
Max. capacity: 24 × 1.5/2.0 mL tubes
- **Rotor for 15 mL/50 mL conical tubes:**
Max. speed: 7,500 × g (7,800 rpm)
Max. capacity: 6 × 15/50 mL tubes

Real Time PCR Step One Plus (Applied Biosystems, USA)

Applied Step One plus is a 4-color, 96-well PCR instrument; allowing the possible use of PCR strips or 96-well plates for high number of samples.

The System supports various real-time PCR applications such as:

- Gene expression profiling
- SNP genotyping analysis
- Non-coding RNA (e.g. MicroRNA and LncRNA) expression analysis
- The system provides high-quality qPCR results



Triple Quadrupole LC-MS/MS (Shimadzu, Japan)

- The triple quadrupole LC-MS/MS equipped with a heated Electrospray ionization (ESI) probe, offers superior performance and sensitivity to meet the demands of quantitative analysis in complex matrices including biological fluids, food, and environmental samples.
- The triple quadrupole mass analyzer can be used for structural elucidation, as it provides information about fragmentation patterns.
- The system is equipped also with a photodiode array (PDA) detector for RP-HPLC/PDA analysis.